

Exp. No. (1)

Preparation of Poly Methyl Methacrylate (PMMA) by Addition Polymerization

Introduction


A polymer is a long-chain molecule that is composed of a large number of repeating units of identical structure.

The repeat unit of polymer is usually equivalent or nearly equivalent to the monomer such as the repeat unit of polyethylene ($-\text{CH}_2-\text{CH}_2-$)_n is $n\text{CH}_2=\text{CH}_2$.

The number of repeat unit along the polymer chain is called the degree of polymerization.

Monomer is the compound or compounds, used in the preparation of polymer.

For example ethylene is a monomer that can be polymerized to polyethylene:



Two monomer molecules first react together to form a dimer. The dimer may then react with a third monomer to yield a trimer and so on.


It is convenient to classify polymers according to the kinetics of the polymerization:

1. Chain-growth polymerization (addition polymerization).
2. Step-growth polymerization (condensation polymerization).

31
Chain-growth thpolymerization is that attachment of a monomer to an (active) chain. The active end may be a free radical or an ionic site (i.e., anion or cation). High-molecular-weight polymer is formed in the early stages of a chain-growth polymerization.

Chain-growth polymerization has three principal steps:



1. Initiation process, in which active monomer is formed




2. Propagation or growth of the active chain by sequential addition of monomers




3. Termination of the active chain to give the final polymer product by different reaction



Initiator is usually molecule that can be decomposed thermally or by irradiation to yield a pair of initiating radicals, such as peroxides, azo compounds and redox system.





Azobisisobutyronitrile (AIBN)

Redox system



Apparatus

1. Balance
2. Heavy-walled glass tube (test tube)
3. Beaker
4. Watch glass
5. Thermometer
6. Oven

Reagent and materials

1. Methyl methacrylate MMA (monomer)
2. Benzoyl peroxide (Initiator)
3. Ethanol (precipitator)

Procedure

1. Put 2 ml of methyl methacrylate individually in two tubes.
2. 0.03g of benzoyl peroxide is added to each tube and shake gently, the tubes are heated inside water bath at 80°C for 15 min.
3. When the contents of the tubes become viscous due to polymerization process, one of MMA tubes is transfer to dry beaker.
4. 10ml ethanol is added to the content of the beaker gradually with shaking, and the precipitated is collected and transferred to the watch glass then dried in oven at 40°C and weighted .
5. The second tube of MMA monomer is heated for 20min, and the same procedure of polymer precipitation in first tube is applied to the second tube and the precipitate is collected, dried and weighted.
6. Calculate the yield percentage of the PMMA prepared by the two procedures.

Exp. No. (2)


Preparation of Polyester by Condensation Polymerization

Introduction


Step-growth polymerization is a random reaction of two molecules that may be any combination of a monomer, oligomer, or a longer-chain molecule, and usually small molecules like H_2O , CH_3OH , HCl , etc. are eliminated. High-molecular-weight polymer is formed in the early stages


An example of this type of polymerization is the formation of laminating resin polyester (Glyptal) from the reaction of glycerol and phthalic acid.

Polyesters are being commercially produced from the combination of two different monomers, one have two (or more) hydroxyl groups and the other have two (or more) carboxyl groups:




The reaction between glycerol and phthalic anhydride is basically esterification, and H_2O molecule is eliminated.







Phthalic anhydride

+




Glycerol



Linear poly ester

230



Cross linked poly ester

Apparatus

1. Pyrex test tube
2. Thermometer
3. Crucibles

Reagent and materials

1. Phthalic anhydride
2. Glycerol
3. Conc. Sodium hydroxide solution
4. Conc. Sulfuric acid
5. Solvent (chloroform, acetone and ethanol)

Procedure

1. 2.5g of phthalic anhydride and 1g of glycerol are mixed in dry Pyrex tube with stopper and thermometer.
2. Heat slowly with careful stirring mixture inside the hood.
3. Rise the temperature to about 160°C and keep the reaction for 15 min at this temperature until all the solid is dissolved.
4. Some drops of produced polymer is pull in two crucibles for analysis.
5. The temperature of polymerization reaction is then rised slowly until it become 200-240°C where cross-linked polymer is formed.
6. Another few drops of product cross-linked polymer is pull in another two crucibles for analysis.
7. The collected drops of polymer samples, two crucibles contain linear and two crucibles contain cross-linked polyester are tested as following:

a) *Conc. Sodium hydroxide test:*

3ml of conc. Sodium hydroxide is added to one linear polyester crucible and another 3ml for one cross-linked polyester crucible.

b) *Acid test:*

3ml of conc. H_2SO_4 is added to one linear polyester crucible and another 3ml for one cross-linked polyester crucible.

c) *Melt test:*

Heating the prepared polymer samples (linear and cross-linked polyester) on a heater and recorder your notes.

d) *Solubility test:*

3ml of one of the following organic solvents as, chloroform, acetone, CCl_4 , ethanol, or any available organic solvent in the lab.

8. Record your notes and obtained results.

* الراتينجول (تقوم بالقيام بحل البكرة و الترسيت
 * وسط البكرة يكون حبيبات خضراء حل البكرة
 * يكون التفاعل باعث للحرارة (منه الهجوم التوسيعي دالة الحرارة)
 * التزوج الزائدة للتخليط
 * يفضل دالة لزوجة عالية ليسيت عدم استخدام المذيب

Exp. No. (3)

Preparation of Synthetic Rubber (Thiokol)

* معونات البكرة
 * هجوم البكرة على الحوائط والقلاص فيها
 * الزائدة
 * هجوم البكرة على التراكيب لانه لا يتلصق حبات

Introduction

Thiokol rubber is a kind of rubber that is resist the natural oxidation by oxygen, ozone and organic solvents such as oils and gasoline. This makes it useful for engine O-rings, gaskets and hoses which may come in contact with oils.

The preparation of Thiokol rubber involve two steps process:

The first step is the preparation of sodium polysulfide by the reaction of sulfur (S_8) with the strong base (NaOH)



At room condition, sulfur is normally in the form of S_8 rings and chain.

The reaction with sodium hydroxide produces a mixture of mostly Na_2S_8 with other chain lengths of sulfur presents, this mixture is called polysulfide.

The second step of preparation is the reaction of sodium polysulfide with 1,2-dichloroethane.



Apparatus

1. Conical flask
2. Graduated cylinder
3. Beaker
4. Water bath

Reagent and materials

1. Sodium hydroxide solution (20%)
2. Powdered sulphur
3. 1,2-dichloroethane
4. Distilled water

Safety

1,2-dichloroethane is flammable, avoid any fires or sparks. It is also an irritant to the skin and eyes and it is toxic by ingestion, work only under a fume hood.

Disposal

Dispose all liquid waters in the proper waste bottle provided.

Procedure

1. 3g powder sulphur is keep in clean and dry round bottom flask.
2. 12ml (20%) NaOH solution is added to the sulfur powder.
3. 40ml distilled water is added to the mixture
4. Reflux the mixture for 30 minutes.


- صافي
5. Keep the reflux contents to cool, filter the remaining mixture to remove the unreacted sulfur.
 6. 15ml of the red brown filtrate solution of Na_2S_x is added to another round bottom flask contain 5ml of dichloroethane.
 7. Heat the mixture at 70°C use water bath with stirring for 45 minutes, the synthetic yellow rubber (Thiokol) is formed.
 8. Cool the contents and wash the product well with distilled water.
 9. Dry and weight the product.

Exp. No. (4)

Preparation of Poly (vinyl alcohol) (PVA)

Introduction

The monomer vinyl alcohol ($\text{H}_2\text{C}=\text{CH}-\text{OH}$) is unstable and it transforms into the tautomeric acetaldehyde and it must be prepared by indirect methods, however the poly (vinyl alcohol) is prepared by the alcoholysis of poly (vinyl alcohol) (the less accurate term "hydrolysis" and "saponification").



Ethanol or methanol can be used to effect the alcoholysis with either acid or base as catalyst.

Apparatus

1. Round bottom flask
2. Glass condenser
3. Water bath
4. Buchner funnel

Reagent and materials

1. Poly (vinyl acetate)
2. Methanol
3. Potassium hydroxide
4. Distilled water

Procedure

1. Set up reaction flask on a water bath.
2. Charge it with 0.6g of poly (vinyl acetate) (PVAc).
3. 15ml of methanol added to the flask.
4. Fit the flask with a reflux water condenser attached to a CaCl_2 drying guard tube.
5. Reflux the contents until PVAc dissolved (25-30) minutes.
6. Add 5ml of methanol containing 0.4g of potassium hydroxide when the contents of the flask are still warm.
7. Reflux the mixture further for 1 hour.
8. Poly (vinyl alcohol) formed due to the hydrolysis of PVAc.
9. Collect the product using Buchner funnel and wash it with 5ml of methanol.
10. Dry the product and weight it.
11. Report the yield.

Exp. No. (5)

Preparation of Urea-Formaldehyde Resins

Introduction




Urea-formaldehyde resins is prepared by step-growth polymerization (condensation polymerization).

The resin is used as adhesive for the bounding of plywood and other structure wood products, because urea formaldehyde resins are lighter in color than phenol formaldehyde resins, they are used for interior plywood and decorative paneling.

The reaction between urea and formaldehyde involved two steps, The first step is the addition of formaldehyde to urea to form methylolated urea, the second step is the condensation between methylol and amino group to form methylene-linked polymers or between two methylol groups to form cross-linked polymers.

As in phenolic resins, the formation of urea formaldehyde resins is known to be greatly influenced by both acid and base catalysts.

The condensation of methylol unites to yields polymer is shown in the following reaction:



Apparatus

1. Beaker
2. Round bottom flask
3. Condenser
4. Buckner funnel
5. Water bath
6. Thermometer

Reagent and materials

1. Urea
2. Formaldehyde
3. NaOH (dilute)
4. HCl (dilute)
5. Filter paper
6. Litmus paper

Procedure

1. Weight 3.5g of urea and crushed it in dry beaker.
2. Add 10ml of formaldehyde solution and stirr until solid urea is dissolve.
3. Few drops (4-5 drops) of sodium hydroxide solution is added for solution medium to become basic, use litmus paper as indicator.
4. Pour the mixture into dry and weight round bottom flask with coke fit condenser.
5. Heat the mixture in water bath for 1 hour at 70-80°C use thermometer to keep the temperature constant.
6. Cool the mixture and add few drops (4-5 drops) of dilute hydrochloric acid to neutralize the mixture (pH = 6-7). Use litmus paper as indicator.
7. Filter the solution using Buckner funnel. Dry the product and weight it for yield percentage.

المادة الموجودة في التذيق التي ماتحتوي السكر للحيات (حقلون) وتختلفها بالقيمة
بروي الكيفيات

* لا يذوب في الماء والاكحول ولا في ليزر في 100% و 112% واحدة فيزيك له صيغة $Cu(H_2O)_3(OH)_2$


Exp. No. (6)

Preparation of Rayon (Pure Cellulose)

Introduction

Cellulose is among the simplest of the natural polymers, consist of

D-glucose repeat unit linked through carbon 1 and 4 by β -linkages.



Cellulose هي بوليمر طبيعي

الغذاء الرئيسي من حيث مصادر البوليور هي مادة المواد الطبيعية والبكرتية الموجودة في القطن الخشام

Cotton and wood is the most important source of cellulose. Wood

has now become the main source of cellulose, the substance must be extracted by pulping.

الخشب والورق
السليلوز هو
المستخلص
من الخشب والورق

هناك ثلاثة أنواع من السليلوز

There are three kinds of cellulose; α -, β -, γ -cellulose:

1. α -cellulose: it is the part which is not dissolved in 18% NaOH solution, it is called pure cellulose, which has high molecular weight and long chain. *سلسلة طويلة غير متفرقة*
2. β -cellulose: it is part which is dissolved in 18% NaOH solution, it can be precipitated by acid and has low molecular weight.
3. γ -cellulose: it is dissolved in 18% NaOH solution and cannot precipitated by acid.

Cellulose can be prepared in two processes:

1. Copper ammonium process.

2. Viscous process.





a- Steeping process
b- Viscous Preparation
c- Pure cellulose fiber preparation

The viscous process cellulose is soaked in caustic soda solution for about 3hr. the alkaline solution is removed and the product is treated with CS_2 this gives cellulose xanthate, which is dissolved in $NaOH$ solution to give viscous solution. The solution filtered and forced through a spinner into a dilute H_2SO_4 solution.

In this experiment pure cellulose can be prepared by viscous which includes:

- Steeping process.
- Viscous preparation.
- Pure cellulose fiber preparation.

Reactions



Note: Cellophane production

When viscose is extruded through a slit die into an acid bath, a film of regenerated cellulose known as *cellophane* were results.

Apparatus

1. Beaker (250ml)
2. Glass rod
3. Cock fit conical with stopper
4. Syringe

Reagent and materials

1. Filter paper or cotton
2. NaOH solution
3. Carbon disulfide CS_2
4. Dilute H_2SO_4
5. Sodium bisulfite Na_2SO_3
6. Zinc sulfate $ZnSO_4$

Procedure A

2. 18% NaOH
1. Put 1g of filter paper or cotton in conical flask or beaker.
 2. 25ml of 18% NaOH solution is added to the conical flask.
 3. Leave the conical flask contents for 1hr use stirring from time to time.
 4. The contents of the conical flask is squeeze between two watch glass, then the squeezed filter paper is cutted into small pieces.

5. Put the pieces in dry conical flask with stopper and left it to the next week.

Procedure B

1. Add 12.5ml of CS_2 to the conical flask which containing cutted pieces of filter paper (prepared last week).
2. Close the conical flask by stopper and keep for 1 hr with stirring from time to time until the color of the filter paper became tallow-red.
3. Pour the excess of CS_2 and squeeze the pieces of paper in mortar and 3-5 drops of 4% of NaOH solution is added to prepare viscous solution with yellow-red color.

Procedure C


1. Spin bath is prepared by mixing the following materials.
 - a) 5ml of H_2SO_4
 - b) 6.2g of Na_2SO_4
 - c) 0.5g of ZnSO_4
 - d) 34.5ml of water
2. Inject the viscous solution into the spin bath to produce regenerated cellulose fibers.
3. Remain the cellulose fibers in spin for 10sec.
4. The prepared rayon filaments are washed, dried and weighted.

Exp. No. (7)

Preparation of Poly (Hexamethylene Sebacamide), Nylon 610 by Interfacial Polycondensation Method

Introduction

Nylon 610, poly (hexamethylene sebacamide), is formed by a simple method and at room temperature called interfacial poly condensation method.



In interfacial polymerization, two complementing monomers (for example, a diamine and a diacid chloride as in this experiment) are dissolved in two immiscible solvents, neither of which is a solvent for the resulting polymer. When the solution are mixed, a film of high polymer is immediately formed at the interface. This can be pulled out as a continues fiber (the "nylon rope trick"); or, if the solution are mixed with violent agitation to produce a large surface area, polymer can be formed in high yield very rapidly. This nonequilibrium method eliminates three of the disadvantages of bulk polymerization: the need for high temperatures, long reaction times, and exact stoichiometric equivalence.

21

Warning

The reaction need steam of nitrogen to prevent explosion may occur in presence of oxygen.

Apparatus

1. Separator funnel (300ml).
2. Large fritted-glass filter, 45-inch diameter, or equivalent.

Reagent and materials

1. Hexamethylene diamine, reagent grade.
2. Sebacoyl chloride
3. Sodium hydroxide
4. Carbon tetrachloride (CCl_4)
5. Acetone
6. Methanol

Procedure

1. Dissolved 2.32g hexamethylene diamine and 1.6g sodium hydroxide in 330ml distilled water, and place it in blender.
2. Dissolve 4.78g sebacoyl chloride in 250 ml CCl_4 and place in separator funne;. Support the separator funnel with its stem passing through the cover of the blender (beaker have magnetic bar).
3. Turn the blender on to maximum speed, and add the contents of the separator funnel rapidly (15 sec.). continue to stir vigorously for about 3 min.
4. Stop the blender, and collect the polymer on the fritted-glass filter.

5. Return the polymer to the blender, and wash with 500ml of 1:1 water : methanol for 2-5 min. then repeat step 4.
6. Repeat step 5 using a water : acetone mixture.
7. Dry the polymer in a vacuum oven for 24 hr. at about 80°C.
8. The produced polymer could be used for primary identification of polymers.

Exp. No. (8)

Preparation of Phenol-Formaldehyde resin (Bakelite)


Introduction

Phenolic resin was the first totally synthetic coating material. In 1909, the Bakelite was the first phenolic resin produced commercially and it was replace for cellulose nitrate.

The base-catalyzed addition of formaldehyde to phenol would form *Resole formation* and acid-catalyzed reaction would form *Novolac formation*.

In resole formation the reaction of phenol and excess formaldehyde under basic conditions would form mono-, di-, an trimethylolphenols.

The resole formation can also be obtained by an initial acid-catalyzed reaction of formaldehyde with excess phenol.



Apparatus

1. Beakers
2. Conical flask
3. Glass rod
4. Water bath
5. Cylinder

Reagent and materials

1. Glacial acetic acid.
2. Formaldehyde solution 40%
3. Phenol
4. Concentration hydrochloric acid
5. Distill water

Procedure

1. Pour 5ml of glacial acetic acid in a beaker.
2. Add 2.3ml of formaldehyde solution in same beaker.
3. Add 2g of phenol with stirring and heating (phenol should be stirred in the formaldehyde solution for 10 minutes to make saturated solution).
4. Add 2-3ml of concentrated HCl to the mixture (concentrated HCl should be added stepwise with continuous stirring for heat dispersion).
5. Filter the large piece of pink resin from its mixture.
6. Wash several time with distilled water, dry the resin and calculate its yield percentage.

Precautions

1. The reaction sometimes is vigorous and it may jump from the beaker, be aware the beaker while adding conc. HCl.
2. Phenol should be stirred with formaldehyde solution for 10min. for saturated solution formation.
3. Conc. Hydrochloric acid should be added drop wise carefully.

Exp. No. (9)

Successive Precipitation of Polymers

Introduction

Most of polymers are consisted of different chain length. So there are differences in their molecular weights. The polymer in this case is called (polydispersed) polymer, while polymer which have chain with similar in their molecular weight is called (monodispersed) polymer. Separation of mixture of molecules with different or similar molecular weights are impossible by the ordinary methods, like distillation or crystallization. Separation is often done by methods depend upon the ability of polymer chains to dissolve according to their molecular weights (their ability to dissolve is decrease as their molecular weight is increase). The separated parts with almostly similar range of molecular weight are measured in their molecular weights. The successive methods followed for separation of polymers according to their chain lengths:

1. Precipitation separation
2. Separation by using different solvent
3. Gradual thermal
4. Column separation

The method used in this experiment is precipitation method which depends on precipitating agent added to the polymer solution while the polymer molecules are not homogenous in their molecular weights so it precipitate gradually according to their molecular weights (this mean big molecular precipitate first then it is followed by the smaller one and so on).

Apparatus

1. Burette 25ml
2. Beaker 150ml
3. Glass rod

Reagent and materials

1. Polystyrene
2. Butyl acetate
3. Ethanol or methanol

Procedure

1. Dissolved 2.5g of polystyrene in 25ml *n*-butyl acetate.
2. Fill burette with precipitating agent like ethanol or methanol.
3. The precipitating agent is added stepwise from the burette with continuous stirring till white precipitate is appeared.
4. Stop the addition of the precipitating agent and separate the white precipitate.
5. The precipitating agent is further added as drops to the separated white precipitate for extra collection.
6. The remained clear solution is treated further with precipitating agent for collection of second part of polymer precipitate, and the procedure is repeated till three different polymer chain lengths are collected.

Each part is putted alone in individual crucible and left next week to dry, then its molecular weight will be calculated.


Exp. No. (10)

Viscosity Average Molecular Weight Determination

Introduction

Measurement of the viscosity of dilute polymer solution provide the simplest of most widely used technique for routinely determining molecular weights for polymer samples.

However, it is not an absolute method because each polymer system must first be calibrated with absolute molecular weight determinations (usually by light scattering) run on fractionated polymer samples. Viscosities are measured at concentrations of about 2 g dL^{-1} of solvent by determining the flow time of a certain volume of solution through capillary of fixed length. Flow time in seconds is record as the time for the meniscus to pass between two designated marks on the viscometer as shown in fig (1). Viscosities of the polymer solutions are measured at constant temperature usually at $30.0 \pm 0.1^\circ\text{C}$.



Ubbelohde viscometer

Thus concentration can be reduced without having to empty and refill the viscometer, the additional advantage is that dilution of solution can be done in the viscometer itself.

Caution: It is necessary to ensure that the polymer solution is free from any solid particle, either due to dust or incompletely dissolved polymer, as any particle that got stuck at capillary tube would affect the flow time.

Viscosity can be expressed in several ways as described below.

Relative viscosity (η_r)

Relative viscosity is the ratio of solution viscosity to solvent viscosity, which is proportional to a first approximation for dilute solutions to the ratio of the solution flow time to solvent flow time.

$$\eta_r = \frac{\eta}{\eta_0} = \frac{t}{t_0} \quad (1)$$

Where η and η_0 are the viscosities of the solution and solvent respectively, t and t_0 are the respective flow time of the solution and the solvent respectively.

Specific viscosity (η_{sp})

$$\eta_{sp} = \frac{\eta - \eta_0}{\eta_0} = \frac{t - t_0}{t_0} = \eta_r - 1 \quad (2)$$

Both η_r and η_{sp} are dimensionless (without a unit). As concentration increases, so does viscosity. Hence to account for the concentration

effects, the specific viscosity is divided by concentration to give the reduced viscosity (η_{red})

$$\eta_{red} = \frac{\eta_{sp}}{C} \quad (3)$$

Where concentration C is commonly expressed as g per 100 ml.

Apparatus

1. Ubbelohde viscometer
2. Water bath
3. Thermometer
4. Beaker
5. Pipet (10ml)
6. Stop watch

Reagent & material

1. Polystyrene
2. Benzene

Procedure

(a) Determination of t_0

1. Set the water bath at $30.0 \pm 0.1^\circ\text{C}$. Measure 10 ml of benzene (which has been filtered through a filter stick) into an ubbelohde viscometer. Allow the system to attain thermal equilibrium (let the viscometer be in the water bath for at least 2 minutes).

2. Allow the solution to flow under gravity. Start timing with a stopwatch when the upper meniscus reaches A, and stop the timing when the meniscus reaches B. Repeat 2 to 3 times, and take the average of the three most consist readings to be the solvent flow time t_0 .
3. Pour out all benzene and hang the viscometer upside down in the fume cupboard to drip dry while you prepare the polystyrene solution as follow.
4. Prepare a solution of the polystyrene sample (0.2g/10ml) by weighing accurately 0.2 g of polystyrene into the beaker. Add 10ml of benzene into the beaker and allow the polystyrene to dissolve completely (need to wait at least 15 minutes, with occasional shaking).

(B) Determination of solution flow time t

5. Pour the prepared polystyrene solution into the viscometer. Follow the same procedure as in (a) to measure the solution flow time (t).
6. Dilute the concentration of the solution in the viscometer by adding 2 ml of solvent directly into the viscometer. Make sure you shake the viscometer for a few minutes to ensure homogeneous mixture, before you start to determine the solvent flow time.
7. Repeat the dilution 2 to 3 times, each time by adding 2ml solvent into the viscometer. For each concentration, repeat measurement of flow time at least 3 times and the average value is computed.

Note: At the end of the experiment, pour the content of the viscometer into a waste bottle. Fill the viscometer with clean solvent, before you return it back to the laboratory assistance.

Analysis of data

Tabulate your results as follow:

Measurement of solvent flow time t_0

No.	Solvent flow time t_0/s
1	
2	
3	


Average of three most constant measurement

gdL ⁻¹	Solution flow time (sec)			Average t (sec)	$\eta_r = t/t_0$ Relative viscosity	$\frac{t}{t_0} - 1 = \eta_{sp}$ Specific viscosity	η_{sp}/C Reduced viscosity
	1	2	3				

Plot (η_{sp}/C) versus (C) in graph paper, and determine the intrinsic viscosity $[\eta]$. Determine the viscosity average molecular weight of your polystyrene sample.

The intrinsic viscosity of a polymer is related to the viscosity average molecular weight by the Mark-Houwink Equation.

$$[\eta] = K (M_v)^a$$



Where K and a are constant for a specific polymer and are dependent on the temperature and solvent.

Table 1: The viscosities constants for some important polymers


Polymer	Solvent	Temp (°C)	$K \times 10^{-4}$	a
Polystyrene	Benzene	20	1.23	0.72
	Benzene	25	4.17	0.60
	Benzene	30	1.06	0.735
	Toluene	20	0.007	0.93
	Toluene	25	1.17	0.72
	Toluene	30	1.1	0.725
Poly Methyl methacrylate	Acetone	30	14.7	0.52
	Benzene	25	0.64	0.76
	Toluene	25	0.71	0.73
	<i>n</i> -butylacetate	25	2.5	0.68

Exp. No. (11)**Determination the Degree of Crystallinity****Introduction**

Crystallization of polymers have a ^{کثیر}significant effects on optical, mechanical, thermal and chemical properties of polymers.

The degree of crystallinity is estimated by different analytical methods and it typically ranges between 10 and 80%, thus crystallized polymer are often called (semicrystalline).

The properties of semicrystalline polymers are determined not only by the degree of crystallinity, but also by the size and orientation of molecular chain.

**Amorphous****Semi-Crystalline**

The arrangement of molecular chains in amorphous and semicrystalline polymers

Most methods of evaluating the degree of crystallinity assume a mixture of perfect crystalline and totally disordered areas; the transition

areas are expected to amount to several percent. These methods include *density* measurement, *differential scanning calorimetry* (DSC), *X-ray diffraction* (XRD), *infrared spectroscopy* and *nuclear magnetic resonance* (NMR). The measured value depends on the method used, which is therefore quoted together with the degree of crystallinity.

Apparatus


1. Volumetric flask (5ml).
2. Glass rod.
3. Graduated cylinder.

Reagent and materials

Many samples of crystalline and amorphous polymers such as PE, PS, PP, PET and PVA.

Procedure

1. Put 5ml of distill water in cylinder 10ml capacity.
2. Put small piece of unknown polymer in the cylinder.
3. Notice the piece of polymer floated on the water surface or sink:



4. The quantity of solution is weighted with volumetric flask (W₂).
5. Calculate the density of solution which equal the density of polymer as

$$\text{Density} = \frac{W_2 - W_1}{V}$$

Where W₁ = weight of volumetric flask and $\frac{V}{2}$ = volume of volumetric flask.

6. Determined the degree of crystallinity of polymer by

$$\% \text{ Crystallinity} = \frac{P - P_{am}}{P_{cry} - P_{am}} \times 100$$

Where P = density of polymer sample, P_{am} = density of amorphous polymer and P_{cry} = density of crystalline polymer.

7. Discussion of the results.

Table 2: Densities of crystalline and amorphous polymers


Polymer	Crystalline P_{cry}	Amorphous P_{am}	$P_{cry} - P_{am}$
PE	1.00	0.852	0.148
PP	0.937	0.854	0.083
PS	1.11	1.054	0.057
PET	1.455	1.335	0.120
P.V.al	1.345	1.269	0.076

Exp. No. (12)

Techniques of Polymerization

Introduction

Phase systems in polymerization can be classified as:



In this experiment we shall discuss bulk polymerization and suspension polymerization:

The characterization of these methods shown in table 3.

Method	Homogenous
Bulk polymerization	The reaction is carried out within the monomer which is catalyzed by initiator, highly exothermic process, increasing in viscosity. Ex. Preparation of PMMA, PVC, PVA.
	Both the monomer and the polymer produced are soluble in a suitable solvent, easy dissipation of heat evolved as a result of exothermic reaction. The solvent is occluded and firmly traps the polymer for this reason solution polymerization is mainly applied when solutions of polymer are required for technical application for e.g. as lacquers, adhesives. The products of this method are relatively low molecular weight because of the possibility of chain transfer to solvent.
Solution polymerization	Heterogeneous
	The monomer containing initiator, modifier etc. is dispersed in a solvent (water) by vigorous stirring. For a stable suspension of monomer in the solvent some stabilizers are added. These may be protective colloids.

Apparatus

1. Beaker
2. Test tube (3)
3. Thermometer
4. Glass rod
5. Watch glass
6. Water bath

Reagent and materials for

Procedure A

1. Methyl methacrylate monomer
2. Benzoyl peroxide
3. Ethanol

Procedure B

1. Methyl methacrylate monomer
2. Benzoyl peroxide
3. Benzene
4. Ethanol

Procedure C

1. Methyl methacrylate monomer
2. Benzoyl peroxide
3. Distilled water
4. Ethanol

Procedure A

(Bulk polymerization of methyl methacrylate)

1. Add 2.5ml of methyl methacrylate monomer in dry test tube with 0.03g of benzoyl peroxide.
2. Heat the mixture in water bath at 80°C for a time (15min) with stirring.
3. Take the mixture out of water bath, then add (5ml) of ethanol with stirring to precipitate the polymer.
4. Separate the precipitate polymer from the solution, then carry it to weight filter paper.
5. Dry, weight the polymer and calculate the yield percentage.

Procedure B

(Solution polymerization of methyl methacrylate)

1. Add 2.5ml of methyl methacrylate monomer in dry test tube with 0.03g of benzoyl peroxide and 2.5ml of benzene.
2. Heat the mixture in water bath at 80°C for a time (15min) with stirring.
3. Take the mixture out of water bath, then add (5ml) of ethanol with stirring to precipitate the polymer.
4. Separate the precipitate polymer from the solution, then carry it to weight filter paper.
5. Dry, weight the polymer and calculate the yield percentage.

Procedure C

(Suspension polymerization)

1. Add 2.5ml of methyl methacrylate monomer in dry test tube with 0.03g of benzoyl peroxide and 2.5ml of distilled water, 0.1g PVA or CMC.
2. Heat the mixture in water bath at 80°C for a time (15min) with stirring.
3. *Take the mixture out of water bath, then add (5ml) of ethanol with stirring to precipitate the polymer.*
4. Separate the precipitate polymer from the solution, then carry it to weight filter paper.
5. Dry, weight the polymer and calculate the yield percentage.

Note: Add PVA or CMC for protective colloids.

Exp. No. (13)

Primary Identification of Polymers

Introduction


The reaction of various tests are directly related to the structure of the polymeric materials.

A useful link can be made in this way between structure and behaviors.


It should be realized that these simple tests is not always possible to give a positive identification of every type of plastic material. As some plastics are modified by additions, fillers, colors ... etc.

1st test: Thermosetting or Thermoplastic

The material should be cutted to a small pieces, exposing to heat and notice if



Another tests are followed.



2nd test: Density measurement (Specific gravity)

Put small piece of plastic material in cylinder containing water and the floatation.




Table 4: Density and specific gravity of some polymer

	Polymer	Density	Specific gravity
Float on water	PP	0.850	0.90
	LDPE	0.90	0.92
	HDPE	0.96	0.96
Sink in water	PS	1.04	1.07
	ABS ✕	1.10	1.10
	Nylon ✕	1.10	1.10
Sink in hypo	PMMA	1.20	1.20
	PVC	1.40	1.40
	UREA	1.50	1.50
	PTFE	2.00	2.75
	Cellulous	1.28	1.30

3rd test: *Burning test*

The plastic materials can be classified into three types according to their burning characteristic.

1. Burned and did not extinguish after the fire source was keep back (extinguishable) e.g. PS, PE, PP, ...
2. Burned and self extinguish after the fire source was keep back (extinguishable) e.g. PVC, Rubber.
3. Unburned, PTFE.

Procedure

Hold a small in flame. If it ignites, withdraw it from the flame a distance of about 6cm. Make observations as listed below:

1. Degree of flammability: burns; self extinguishing or not, dose not burn.
2. Color of flame: yellow and orange, blue, blue with yellow edge, etc... (notice the tip and the base of the flame).
3. Color of smoke: no smoke; black, white.
4. Type of flame: steady, tall, short, sputtering.
5. Vapor after extinguishing: acidic or basic by litmus paper.

Compare your results with table (5) to decide the type of the polymer.

Table 5: Typical results of burning tests

Polymer	Color of flame	Kind of smoke	Vapor	Other features
PE	Blue with yellow tip	Absent	Neutral	Flaming plastic
PP	Yellow with blue base	Absent	Neutral	Flaming plastic
PTFE	Will not burn	Non	Non	Non
PS	Orange-yellow with blue base	Heavy black smoke	Neutral	Spurting burns
PVC	Yellow with green base	White smoke	Acid	Difficult to burn and self extinguish
Cellulose acetate	Yellow and turn brown at base	Light smoke	Acid	Burn, strong smell of vinegar
PMMA	Blue and yellow bubble at base	Little black smoke	Neutral	Burns and not extinguish
Nylon	Blue with yellow tip	Non	Neutral	Burns and flame lead to go out
Acetal resin	Pale blue	Non	Neutral	Burn
Poly carbonat	Yellow	Black smoke	Neutral	difficult to light
Phenol resin	Yellow with spark	Non	Neutral	Difficult to burn and self extinguish
Urea resin	Yellow with light blue edge	Non	Basic	Burn and flame lead to go out
Melamine resin	Pale yellow with light blue-green edge	non	Basic	Difficult to burn
Poly ester	Yellow with blue base	very dense smoke	Neutral	burn
Poly urethan	Yellow with blue base	Non	Neutral	burn
ABS	Yellow with blue base	Black smoke	Neutral	Difficult to burn and self extinguish

4th test: Melting and odor characteristics

Heat a metal plate over a Bunsen flame or electric plate and piece of polymeric material on it.

Observe the melt behavior (softness, hard, flows, bubbles, swells does not char but melts and chars). For odor, use a description by relating it to a familiar smell (burning rubber, sweet fruity, marigold, burning paperacid, rancid rubber). In most cases a strong odor will be noticed when the material is burning. Check your result against table 6.

Table 6: Melting and odor characteristics

Polymer	Melting behavior	Odor
PE	Melt and bubbles fiercely around the edge	Waxy burning candle
PP	Becomes soft, tends to form a ball	Waxy smell
PTFE	Does not melt	Non
PS	Melt and bubbles around edges	Marigold (styrene)
ABS	Becomes soft but not fluid	Burning rubber
PVC	Softens and forms a ball, chars at the bottom	Characteristic smell
Cellulose acetate	Melt and emits white smoke boils and chars	Vinegar
PMMA	Boils and bubbles	Fruity, floral smell
Poly carbonat	Bubbles develops acellular structures before decomposing	Non
Nylon	Melt round edge and chars	Burning hair or wool
Phenol resin	Does not melt	Carbolic smell
Urea resin	Does not melt	Fishy smell
Melamine resin	Does not melt	Fishy smell

5th test: Solubility test

The solubility process of polymer is different than the solubility of organic materials (small molecules). In polymer solubility, the polymer swells in primary stage. In some cases there is only swelling but no stability. This case is not solubility so the polymer dissolving needs long time with stirring. In this test need two solvents acetone and *n*-propanol.

Procedure

1. Shave off small amount of plastic materials, place the shaving into a test tube.
2. Add just enough solvent to cover the shavings.
3. Stir with glass rod for a few minutes and note the results.
4. If the plastic does not dissolve, leave for many minutes and stir again.

Check your results against table 7.

Table 7: Typical results at solubility tests

Polymer	Dissolved in acetone	Dissolved in <i>n</i> -propanol
PE	×	×
PP	×	×
not (PTFE)	×	×
PS	○	×
ABS	○	×
PVC	+	×
(Cellulose acetate)	○	+
PMMA	○	×
Nylon	×	×
All resin	×	×

× = not soluble

○ = Soluble

+ = slightly soluble


Exp. No. (14)

Emulsion Polymerization

Introduction

In emulsion polymerization, the monomer is dispersed in water containing an emulsifying agent such as a soap or detergent. Initiator radicals, usually of the redox type, are generated in the aqueous phase and diffuse into soap micelles swollen with monomer molecules. As monomer is used up in the polymerization reaction, more monomer migrates into the micelles to continue the reaction. Termination of polymerization occurs by radical combination when a new radical diffuses into the micelle. Because only one radical is present in the micelle prior to termination, extremely high molecular weights are obtainable, generally too high to be of practical value unless compounds called *chain transfer agents* are added that control the degree of polymerization.

Emulsion polymerization is widely used in industry, and is particularly useful for manufacturing water-based (latex) paints or adhesives.



Apparatus

1. Round bottom flask with 2 neck.
2. Condenser.
3. Glass stirrer.
4. Thermometer.
5. Heating bed.

Reagent and materials

1. Potassium thiosulfate.
2. Potassium persulfate.
3. Acrylo nitrile monomer.
4. Aluminum sulphate.

Procedure

1. Charge the following materials in a 100ml two-necked round bottom flask:
 - 5ml of acrylo nitrile monomer.
 - 0.5g of liquid soap.
 - 0.03g of potassium thiosulfate.
 - 0.03g of potassium persulfate.
 - 30ml distilled water.
2. Connect the instruments as in the following scheme.
3. Heat in bath water at 70°C with vigorous stirring for 1hr.
4. Add 5ml of aluminum sulfate solution (2.5%).
5. Filter and wash with warm water and methanol.
6. Dry at 60°C , weight and calculate the percentage of the yield.

