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

Extraction of Genomic DNA From Whole Blood.

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Ammonium Chloride
2.	Potassium Hydrogen Carbonate
3.	SDS
4.	Di, sodium ethylene diamin tetra acetic acid (Na ₂ EDTA)
5.	Tris-Hcl
6.	Isopropanol
7.	Hydrochloric acid
8.	Ethanol
9.	Sodium chloride
10.	Sodium hydroxide
11.	Proteinase k

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
6.	Biohazard Container
7.	Pipette

The Preparation of Buffers:

1. 10M NaOH :

It is prepared by dissolving **40g** of **NaOH** in 80ml of distill water the volume is made up to **100ml** by distill water.

2. R.B.Cs Lysis buffer :

It is prepared dissolving **4.145 g of Ammonium Chloride** , and **0.5 g of Potassium Hydrogen Carbonate**, and **0.185 g of Na₂EDTA** in **450 ml** of distill water, the **pH** is adjust to (**7.6**) using (**10 N**) **NaOH**. The volume is made up to **500 ml** by distill water and sterilize by autoclaving and store **at 4°C**.

3. 10% SDS Solution :

It is prepared by dissolving **100g** of **SDS** in **900ml** of distill water with heating at **68C** to assist dissolution adjust volume **to one liter**.

4. Cell Lysis Buffer:

It is prepared dissolving **1.576 g of Tris - HCL** , and **0.744 g of Na₂EDTA**, and **2.922 g of NaCl** in **450 ml** of distill water, the **pH** is adjust to (**7.4**) using (**10 N**) **NaOH**. The volume is made up to **500 ml** by distill water and sterilize by autoclaving and store **at 4°C**.

5. Proteinase k (20mg/ml):

It is prepared by dissolving 20 mg of Proteinase k in 1 ml of distill water mixing well and store at -20 °C.

6. 5.3 M NaCl:

It is prepared by dissolving 77.44 g of in 250 ml of distill water mixing well and store at Room Temperature.

7. Ethanol (70%):

comprise of **70 ml** absolute ethanol and **30 ml** of distill water.

8. Isopropanol:

Ready to use.

The Procedure:

1. **300 µL** of the thaw or fresh blood sample , is mixed well with **900 µL** of **2X lyses buffer (Buffer A)** in a sterile Eppendorf tube, by inverting several times and the mixture is left on ice or **0°C** for **15 minutes** inverting occasionally.
2. The mixture is then centrifuge in a bench centrifuge at **10.000 rpm** for **10 minutes** to obtain a nuclear pellet.
3. Add a further **600 µL** of **Buffer A**, mix vigorously and centrifuge at **10.000 rpm** for **5 minutes** (this step is necessary when the cell pellet remains red or brown after the first spin).
4. To the remaining white pellet, **500 µL** of **(Buffer B)** was added, and followed by **50 µL** of **10% SDS**, mixed by vortexing vigorously for 1 min until the mixture was re-suspended and well dissolved by gentle pipeting.
5. To Latter mixture **10 µL** of fresh, refrigerated **Proteinase k solution** (20mg/ml) was added and the tubes mixed well by inverting several times and incubated in a water bath at **55 C°** for 90 min, gently invert the tubes a few times during incubation.
6. Following incubation the tubes were placed on ice to cool for 2-5 minutes, **500 µL** of **5.3 M Nacl** was added and vortex for 1 minutes, followed by centrifugation at **15.000 rpm** for **5 minutes** at room temperature.

7. The supernatant is then removed to another Eppendorf tube and **equal volume** of **cold isopropanol** was added.
8. The Eppendorf tube left at **0°C for 1–5 minutes**.
9. The tubes were inverted **5-6 times** gently to **precipitated DNA**, then centrifuged at **10.000 rpm** for **5 minutes**. The supernatant from this last step discarded and **ethanol (70%) was added**, again tubes centrifuged **10.000 rpm** for **5 minutes**.
10. The supernatant discarded and **200-300 µL** of distil water was added to re-suspended (dissolved) **DNA**.
11. The stock **DNA** is kept frozen at **(-20°C)** until use.

Rapid Extraction of High Quality DNA from Whole Blood

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Vortex

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	sucrose
2.	MgCL ₂
3.	SDS
4.	Triton X-100
5.	Di, sodium ethylene diamine tetra acetic acid (Na ₂ EDTA)
6.	Tris-Hcl
7.	Sodium citrate
8.	Hydrochloric acid
9.	Ethanol
10.	Sodium hydroxide

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Tips with different sizes
2.	Eppendorf tubes (1.5 size)
3.	Racks for Eppendorf tubes
4.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
5.	Biohazard Container
6.	Pipette

Blood collection:

Blood should be collected in EDTA-containing vacutainer tubes. As with all body fluids, blood represents a potential biohazard, thus care should be taken in all steps requiring

handling of blood. If the subject is from a known high-risk category, additional precautions may be required. Blood samples can be stored at room temperature for DNA extraction within the same working day or at refrigerator for later uses.

Standard chemicals:

This method uses standard chemicals that can be obtained from any major supplier; we used chemicals supplied by Sigma Co. as follow:

1. EDTA (0.5 M), pH 8.0:

Add 186.1 gr of anhydrous EDTA to 800 ml of distilled water. Adjust pH to 8.0 with NaOH pellets. Make up to 1 liter with distilled water. Autoclave at 15 p.s.i. for 15 min.

2. 1 M Tris-HCl, pH 7.6:

Dissolve 121.1 gr of Tris base in 800 ml of distilled water. Adjust pH with concentrated HCl. Allow mixture to cool to room temperature before finally correcting pH. Make up to 1 liter with distilled water. Autoclave at 15 p.s.i. for 15 min.

3. Red blood cell lysis buffer:

0.01 M Tris-HCl pH 7.6, 320 mM sucrose, 5 mM MgCl₂, 1% Triton X 100.

Add 10 ml of 1 M Tris, 109.54 gr of sucrose, 1.01 gr MgCl₂, adjust pH to 8.0 and finally add 10 ml of Triton X-100 to 800 ml of distilled water, and make up to 1 liter with distilled water. Autoclave at 15 p.s.i. for 10 min. Sugars at high temperature can cause caramelization (browning), which degrades the sugars [5].

4. Nucleic lysis buffer:

0.01 M Tris-HCl, 11.4 mM sodium citrate, 1 mM EDTA, 1 % sodium dodecyl sulphate (SDS).

Take 10 ml of 1 M Tris-HCl (pH 7.6), 3.75 gr of anhydrous EDTA (pH 8.0), 10 gr SDS, 2.94 gr of sodium citrate, and adjust pH to 8.0. Make up to 1 liter with distilled water. Autoclave 15 min at 15 p.s.i.

5. TE Buffer, pH 8.0:

Take 5 ml of 1 M Tris-HCl, pH 7.6, 2 mL of 0.5 M EDTA, pH 8, and make up to 1 liter with distilled water. Adjust pH to 8.0 and autoclave 15 min at 15. p.s.i.

6. Chloroform prechilled to 4°C.

7. Ethanol (100%) prechilled to -20°C.

Procedure of DNA Extraction:

1. Before starting DNA extraction, liquid blood venogects should be shake gently by rotating blood mixer (vortex).
2. Pour 500 µl of blood into a 1.5 ml eppendorf tube and add 1000 µl of red cell lysis buffer. Shaking microfuge tube gently for 3 min (up to homogenizing), then spin for 5 minutes at 7000 rpm.
3. Discard supernatant and repeat steps 1-3 two or three more times to remove hemoglobin. It is important to breakdown the pellet by vortexing and rinses it well in red blood cell lysis buffer in order to clean the white blood cells from residual of hemoglobin.
4. Placing the tube on tissue paper for few seconds downward. Be careful from cross-contamination between different samples.
5. Add 400 µl of nucleic lysis buffer to eppendorf tube. **Note:** if the pellet formed, you must pipette the pellet up to dissolve it.
6. Place tube in water bath for **15 min at 55°C.**
7. Allow to cool to room temperature.
8. Add 100 µl of saturated 5M NaCl and 600 µl of ice-cold chloroform chloroform to eppendorf tube and mix on a rotating blood mixer at room temperature for **15 min** then spin it for 5 minutes at 7000 rpm.
9. Transfer 400 µl of supernatant to a new 1.5 ml tube.
10. Add 800 µl of cold (-20°C) absolute Ethanol and shake it gently then vortex it. DNA should appear as a mucus-like strand in the solution phase.

11. Spin the microfuge tube for one minute at 12000 rpm to precipitate, then discard supernatant carefully and let tube be completely dried in room temperature (Place Eppendorf tube downward on the tissue paper).
12. Add 50µl of TE to it then vortex; keep eppendorf tube of DNA in 4°C or -20°C for later uses. We routinely use about one µl per PCR reaction without adverse affects. DNA can be quantified and diluted to a working concentration at this point or simply use 1 µl per PCR reaction. We expect that the yield of this procedure be 100 to 300 ng/µl, DNA. Using the above method, high quality DNA samples from a sheep population were extracted for gene mapping studies.

Extraction of Genomic DNA From Whole Blood:By Sodium Chlorate:

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The Chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Sodium hydroxide
2.	Sodium Perchlorate
3.	SDS
4.	Di, sodium ethylene diamine tetra acetic acid (Na ₂ EDTA)
5.	Tris-Hcl
6.	Chloroform
7.	Ethanol
8.	Sodium chloride

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20µl, 10-100µl, 100-1000µl)
6.	Biohazard Container
7.	Pipette

The Preparation of Buffers:

1. Red blood cell lysis buffer: 10mM Tris-HCl, 320 mM sucrose, 5 mM MgCl₂, 1% Triton X 100.

It is prepared by dissolving 0.157 g of Tris-HCl, 10.95g of sucrose and 0.1 g of MgCl₂ in 80 ml of distilled water, adjust the pH to 8.0 and 1 ml of Triton X-100 was added, finally the volume is made up to 100 ml by the distilled water and sterilize by autoclave at

15 p.s.i. for 10 min. Sugars at high temperature can cause caramelization (browning), which degrades the sugars.

2. 10 M NaOH :

It is prepared by dissolving **40 g** of **NaOH** in 80ml of distill water the volume is made up to **100 ml** by distill water.

3. Cell Lysis Buffer:

0.4 M Tris-HCl, 150 mM NaCl, 0.06 M EDTA, 10 % sodium dodecyl sulfate (SDS), adjust pH to 8.

It is prepared by dissolving 6.3 g of Tris-HCL and 2.23 g of Na₂EDTA and 0.88 g of NaCl in 80 ml of distill water, the pH is adjust to (8) using (10 N) NaOH. The volume is made up to 100 ml with distill water. Sterilize by autoclaving and add 10 g of SDS.

4. 5 M Sodium Perchlorate:

Dissolve 7 g of sodium perchlorate in 8 mL distilled water, make up to 10 mL.

5. Chloroform prechilled to 4°C.

6. Ethanol (100%) prechilled to -20°C.

7. Ethanol (70%):

comprise of 70 ml absolute ethanol (100%) and 30 ml of distill water.

The Procedure:

1. Place **300 µL** whole blood in a 1.5-mL Eppendorf tube, add **1200 µL** reagent A.
2. Mix on a rolling or rotating blood mixer for **5 min** at **room temperature**.
3. Centrifuge at **10,000 rpm** for **5 min** at room temperature.
4. Discard supernatant without disturbing cell pellet. Remove remaining moisture by inverting the tube and blotting onto tissue paper.
5. Re-suspend pellet in **750 µL** of reagent A by vortexing for **30 seconds** at medium speed. Centrifuge at **10,000 rpm** for **5 minutes**.

6. Add **500 μ L** reagent B and vortex briefly to resuspend the cell pellet.
7. Add **100 μ L** of 5M sodium perchlorate , mix by inverting tube several times.
8. Place tube in water bath for **60min at 65°C**.
9. Allow to cool to room temperature.
10. Add **700 μ L** ice-cold chloroform. Then mix on a rotating mixer for **30 - 60 min**.
11. Centrifuge at **10,000 rpm** for 5 minutes.
12. Transfer upper phase into a clean Eppendorf tube using a sterile pipet.
13. Add equal volume from ice-cold ethanol and invert gently to allow DNA to precipitate. Using a freshly prepared flamed Pasteur pipet, spool the DNA onto the hooked end, transfer to a 1.5-mL Eppendorf tube and allow to air-dry , Or Centrifuge at 10,000 rpm for 5 minutes discard the supernatant.
14. Wash the pellet by add **500 μ L** of **70 %** ethanol. Centrifuge at **10,000 rpm** for **10 minutes**.
15. Resuspend in **200 μ L** TE buffer.

Extraction of Genomic DNA From Animal Tissue:

(Working with Proteinase k)

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Sodium hydroxide
2.	Proteinase k
3.	SDS
4.	Di, sodium ethylene diamin tetra acetic acid (Na ₂ EDTA)
5.	Tris-Hcl
6.	Isopropanol
7.	Hydrochloric acid
8.	Ethanol
9.	Sodium chloride

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Tips with different sizes
2.	Eppendorf tubes (1.5 size)
3.	Glass centrifuge tubes
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20 μ l, 10-100 μ l,100-1000 μ l)
6.	Biohazard Container
7.	Pipette
8.	Petri dishes
9.	

The Preparation of Buffers:

1. 10M NaOH :

It is prepared by dissolving **40g** of **NaOH** in 80ml of distill water the volume is made up to **100ml** by distill water.

2. 10% SDS Solution :

It is prepared by dissolving **100g** of **SDS** in **900ml** of distill water with heating at **68C** to assist dissolution adjust volume to **one liter**.

3. Cell Lysis Buffer:

It is prepared dissolving **1.576 g** of **Tris - HCL** , and **0.744 g** of **Na₂EDTA**, and **2.922 g** of **NaCl** in **450 ml** of distill water, the **pH** is adjust to **(7.4)** using **(10 N) NaOH**. The volume is made up to **500 ml** by distill water and sterilize by autoclaving and store **at 4°C**.

4. Proteinase k (20mg/ml):

It is prepared by dissolving 20 mg of Proteinase k in 1 ml of distill water mixing well and store at -20 °C.

5. 5.3 M NaCl:

It is prepared by dissolving 77.44 g of in 250 ml of distill water mixing well and store at Room Temperature.

6. Ethanol (70%):

comprise of **70 ml** absolute ethanol and **30 ml** of distill water.

7. Isopropanol:

Ready to use.

The Procedure:

1. Place the Petri dish in an ice bucket filled with dry ice, this is a good way to precool the Petri dish before adding the tissue sample.
2. The frozen Animal tissue sample is thawed, put **10-20 mg** of each tissue in a Petri dish containing **100 µL** of cell lysis buffer and chop it into a very small pieces with sterile scalpel blade, then transfer the mixture to sterile Eppendorf tube.
3. To the mixture, **400 µL** of (**Buffer B**) was added, and followed by **50 µL** of **10% SDS**, mixed by vortexing vigorously for 1 min until the mixture was re-suspended and well dissolved by gentle pipeting.
4. To Latter mixture **10 µL** of fresh, refrigerated **Proteinase k solution** (20mg/ml) was added and the tubes mixed well by inverting several times and incubated in a water bath at 55 C° for 90 min, gently invert the tubes a few times during incubation.
5. Following incubation the tubes were placed on ice to cool for 2-5 minutes, **500 µL** of **5.3 M Nacl** was added and vortex for 1 minutes, followed by centrifugation at **14.000 rpm** for **5 minutes** at room temperature.
6. The supernatant is then removed to another Eppendorf tube and **equal volume** of **cold isopropanol** was added.
7. The Eppendorf tube left at **-20°C for 5 minutes**.

8. The tubes were inverted **5-6 times** gently **to precipitated DNA**, then centrifuged at **14.000 rpm** for **3 minutes**. The supernatant from this last step discarded and **ethanol (70%) was added**, again tubes centrifuged **14.000 rpm** for **3 minutes**.
9. The supernatant discarded and **200-300 µL** of distil water was added to re-suspended (dissolved) **DNA**.
10. The stock **DNA** is kept frozen at **(-20°C)** until use.

Extraction of Genomic DNA From Animal Tissue:

(Working with Sodium Chlorate)

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The Chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Sodium hydroxide
2.	Sodium Perchlorate
3.	SDS
4.	Di, sodium ethylene diamin tetra acetic acid (Na ₂ EDTA)
5.	Tris-Hcl
6.	Chloroform
7.	Ethanol
8.	Sodium chloride

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
6.	Biohazard Container
7.	Pipette

The Preparation of Buffers:

1. 10 M NaOH :

It is prepared by dissolving **40 g** of **NaOH** in 80ml of distill water the volume is made up to **100 ml** by distill water.

2. Cell Lysis Buffer:

0.4 M Tris-HCl, 150 mM NaCl, 0.06 M EDTA, 10 % sodium dodecyl sulfate (SDS), adjust pH to 8.

It is prepared by dissolving 6.3 g of Tris-HCL and 2.23 g of Na₂EDTA and 0.88 g of NaCl in 80 ml of distill water, the pH is adjust to (8) using (10 N) NaOH. The volume is made up to 100 ml with distill water. Sterilize by autoclaving and add 10 g of SDS.

3. 5 M Sodium Perchlorate:

Dissolve 7 g of sodium perchlorate in 8 mL distilled water, make up to 10 mL.

5. Chloroform prechilled to 4°C.

6. Ethanol (100%) prechilled to -20°C.

7. Ethanol (70%):

comprise of 70 ml absolute ethanol (100%) and 30 ml of distill water.

The Procedure:

1. The frozen Animal tissue sample is thawed, put **10-20 mg** of each tissue in a Petri dish containing **100 µL** of cell lysis buffer and chop it into a very small pieces with sterile scalpel blade, then transfer the mixture to 2 mL sterile Eppendorf tube.
2. To the mixture, **500 µL** of (**Buffer B**) was added, mixed by vortexing vigorously for 1 min until the mixture was re-suspended and well dissolved.
3. Add **100 µL** of 5M sodium perchlorate , mix by inverting tube several times.
3. The Eppendorf tubes placed in water bath for **90 min** at 65°C, gently invert the tubes a few times during incubation.
4. Allow for 5 min to cool to room temperature.
5. Add **700 µL** ice-cold chloroform. Then mix on a rotating mixer for **30 - 60 min**.
6. Centrifuge at **10,000 rpm** for **5 minutes**.
7. Transfer upper phase into a clean Eppendorf tube using a sterile pipet.
8. Add equal volume from ice-cold ethanol and invert gently to allow DNA to precipitate.
9. Centrifuge at 10,000 rpm for 5 minutes, discard the supernatant.
10. Wash the pellet by add **500 µL** of **70 %** ethanol. Centrifuge at **10,000 rpm** for **10 minutes**.

11. Resuspend in 200 μ L TE buffer.

A simple and rapid leaf genomic DNA extraction method for PCR analysis

Jafar Amani¹, Roohallah Kazemi², Ali Reza Abbasi², Ali Hafez Salmanian¹

IRANIAN JOURNAL of BIOTECHNOLOGY, Vol. 9, No. 1, January 2011

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Tips with different sizes
2.	Eppendorf tubes (1.5 size)
3.	Racks for Eppendorf tubes
4.	Micropipette (transfer pipette) (2-20 μ l, 10-100 μ l,100-1000 μ l)
5.	Biohazard Container
6.	Pipette

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Tris-Hcl
2.	Sodium chloride
3.	Di, sodium ethylene diamin tetra acetic acid (Na ₂ EDTA)
4.	SDS
5.	Chloroform
6.	Isoamyl Alcohol
7.	Isopropanol
8.	Liquid Nitrogen.

The Preparation of Buffers:

1. Extraction Buffer 25 ml:

It is prepared by dissolving **1.0 g** of **CTAB**, **7.87 g** of **Tris-HCL** and **9.3 g** of **Na₂EDTA**, and **14.62 g** of **NaCl** in **40 ml** of distill water. Adjust all to **pH (8)** using **NaOH** and make up volume to **50 ml** with distill water, sterilize by autoclaving.

2. TE buffer: pH 8.0:

Take 5 ml of 1 M Tris-HCl, pH 7.6, 2 mL of 0.5 M EDTA, pH 8, and make up to 1 liter with distilled water. Adjust pH to 8.0 and autoclave 15 min.

3. Chloroform : Isoamyl alcohol (24:1) :

To prepare **100 ml** of this mixture **96ml** of **chloroform** is added to **4ml** of **Isoamyl alcohol** and stored in a dark bottle at **4C** until use.

4. Isopropanol : Ready to use.

The Procedure:

1. In the present study, samples for PCR analysis (usually leaf tissue) were collected into a sterile eppendorf tube. The tissue was softened for 15 s at room temperature, using a sterile tip as grinder, without any buffer.

2. 500 μ l of extraction buffer (200 mM Tris-HCl (pH 7.5), 250 mM NaCl, 25 mM EDTA, 0.5% SDS) (Edwards *et al.*, 1991) was added to the tube and mixed for 5 s.
3. In order to disrupt plant cells completely, the tube was placed in a water bath at 60°C for 60 min.
4. Subsequently, an equal volume of chloroform: iso-amyl alcohol (24:1) was added to the sample, which was then mixed and centrifuged at 15,000 g for 5 min, at 4°C.
5. The aqueous supernatant was transferred to a new eppendorf tube, to which an equal volume of isopropanol was added, mixed and incubated at -20°C for 30 min.
6. Following centrifugation at 15,000 g for 5 min, the pellet was dried and dissolved in (100 μ l) TE buffer.
7. Finally, for precipitation of starch and other insoluble polysaccharides, the tube was placed on ice for 5 min and centrifuged at 15000 g for 2 min; the resulting white pellet was mostly starch (Deshmukh *et al.*, 2007). The supernatant containing the DNA was stable at 4°C.

Plant Genomic DNA Extraction using CTAB:

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	CTAB (Hexadecyl trimethyl-ammonium bromide)
2.	Sodium chloride
3.	PVP (polyvinyl pyrrolidone (vinylpyrrolidone homopolymer)
4.	Di, sodium ethylene diamine tetra acetic acid (Na ₂ EDTA)
5.	Tris-HCl
6.	Chloroform
7.	Ethanol (ice cold).
8.	Sodium chloride
9.	Iso Amyl Alcohol
10.	Liquid Nitrogen.
11.	Ammonium Acetate

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Glass centrifuge tubes
5.	Racks for Eppendorf tubes
6.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
7.	Biohazard Container
8.	Pipette

The Preparation of Buffers:

1. CTAB buffer 50 ml:

It is prepared by dissolving **1.0 g** of **CTAB**, **7.87 g** of **Tris-HCL** and **9.3 g** of **Na₂EDTA**, and **14.62 g** of **NaCl** in **40 ml** of distill water. Adjust all to **pH (8)** using **NaOH** and make up volume to **50 ml** with distill water, sterilize by autoclaving.

5. Ethanol (70%):

comprise of **70 ml** absolute ethanol (100%) and **30 ml** of distill water.

6. Chloroform : Isoamyl alcohol (24:1) :

To prepare **100 ml** of this mixture **96ml** of **chloroform** is added to **4ml** of **Isoamyl alcohol** and stored in a dark bottle at **4C** until use.

7. 7.5M Ammonium acetate :

It is prepared by dissolving **5.782g** Ammonium acetate in **10ml** distilled water sterilize by filtration (micropore filter) with **0.20 µm**.

8. Absolute Ethanol (ice cold): Ready to use.

The Procedure:

1. Grind 100 mg of plant tissue to a fine paste in approximately 500 µl of CTAB buffer.
2. Transfer CTAB/plant extract mixture to a Eppendorf tube.
3. Incubate the CTAB/plant extract mixture for about 10 min at 55°C in a recirculating water bath.
4. After incubation, to each tube add 250 µl of Chloroform : Iso Amyl Alcohol (24:1) and mix the solution on a rotating mixer for 5 min. After mixing, spin the tubes at 14000 rpm for 5 min.
5. Transfer the upper aqueous phase only (contains the DNA) to a clean microfuge tube.
6. To each tube add 50 µl of 7.5 M Ammonium Acetate followed by 500 µl of ice cold absolute ethanol.
7. Invert the tubes slowly several times to precipitate the DNA. Generally the DNA can be seen to precipitate out of solution. Alternatively the tubes can be placed for 10 min at -20 °C after the addition of ethanol to precipitate the DNA.

Following precipitation, the DNA can be pipetted off by slowly rotating/spinning a tip in the cold solution. The precipitated DNA sticks to the pipette and is visible as a clear thick precipitate.

8. To wash the DNA, transfer the precipitate to new microfuge tube containing 500 μ l of ice cold 70 % ethanol and slowly invert the tube. Repeat. ((alternatively the precipitate can be isolated by spinning the tube at 13000 rpm for a 3 min to form a pellet. Remove the supernatant and wash the DNA pellet by adding two changes of ice cold 70 % ethanol)).
9. After the wash, spin the DNA into a pellet by centrifuging at 13000 rpm for 3 min. Remove all the supernatant and allow the DNA pellet to dry (approximately 15 min). Do not allow the DNA to over dry or it will be hard to re-dissolve.
10. Resuspend the DNA in sterile DNase free water (approximately 50-400 μ l H₂O; the amount of water needed to dissolve the DNA can vary, depending on how much is isolated).
11. After resuspension, the DNA is incubated at 65 °C for 20 min to destroy any DNases that may be present and store at 4°C.

DNA Extraction by Modified CTAB Extraction Method According to Moller *et al* :

Global Journal of Molecular Sciences 5 (1): 37-12,2010 India.

ISSN 1990-92-11 IDOSI Publications. 2010

V. Karthikeyan, ²S. Patharajan, ³P. Palani. ¹D. Spadaro, •**ML**. Gullino and ¹A. Garibaldi.

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave

4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Water distiller
8.	Vortex

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Tips with different sizes
2.	Eppendorf tubes (1.5 size)
3.	Glass centrifuge tubes
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
6.	Biohazard Container
7.	Pipette

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Tris-Hcl
2.	Di, sodium ethylene diamine tetra acetic acid (Na ₂ EDTA)
3.	SDS
4.	Proteinase K
5.	CTAB (Hexadecyl trimethyl-ammonium bromide)
6.	Sodium chloride
7.	Chloroform
8.	Isoamyl Alcohol

9.	Ethanol (ice cold).
10.	Isopropanol
11.	sodium acetate
12.	RNase A

The Preparation of Buffers:

1. Lysis Buffer: 10 ml

It is prepared by dissolving **0.156 g** of **Tris-HCL** and **0.04 g** of **Na₂EDTA**, and **200 m g** of **SDS** in **60 ml** of distill water. Adjust all to **pH (8)** using NaOH and make up volume to **10 ml** with distill water, sterilize by autoclaving.

2. Proteinase k (20mg/ml):

It is prepared by dissolving 20 mg of Proteinase k in 1 ml of distill water mixing well and store at -20 °C.

3. 5 M NaCl:

4. It is prepared by dissolving 2.92 g of in 10 ml of distill water mixing well and store at Room Temperature.

5. 10% CTAB Solution :

It is prepared by dissolving **10g** of **CTAB** in **10ml** of distill water with heating at **68°C** to assist dissolution.

6. Chloroform : Isoamyl alcohol (24:1) :

To prepare **100 ml** of this mixture **96ml** of **chloroform** is added to **4ml** of **Isoamyl alcohol** and stored in a dark bottle at **4C** until use.

7. Isopropanol:

Ready to use.

8. Ethanol (70%):

comprise of **70 ml** absolute ethanol (100%) and **30 ml** of distill water.

9. RNase A

It is prepared by dissolving 10mg of RNase A in **1ml** of distill water.

Procedure:

1. 50mg of fungal mycelia was scraped from 10 days old PDA cultures, manually ground in 2 ml of microfuge tubes with micro pestle adding 500 uL of pre-warmed (60°C) TES lysis buffer (100mM Tris pH 8.0; 10mMEDTA pH 8.0; 2% SDS).
2. 5µL (100ug) of proteinase K 20mg/ml were added to the ground material, incubated in 60°C for 60 min.
3. 140uL of 5M NaCl and 64 uL of 10% (w/v) of CTAB were added to the suspension incubated at 65°C for 10 min.
4. DNAs were extracted by adding equal vol. of chloroform: isoamylalcohol (24:1) centrifuged at 14000xg/ 10min.
5. DNA was precipitated by adding 0.6 vol (200 uL)of cold isopropanol and 0.1 vol (100 uL)of 3M sodium acetate pH 5.2 and maintained at -20°C,
6. Centrifuged and washed twice with 70% ethanol suspended in 100ul of TE (10mM Tris pH 8.0; 1mM EDTA pH 8.0).
7. RNA was digested by adding 10mg/ml of RNase A and incubating at 37°C for 45min and stored -20°C for further use.

A Simple and Rapid Method for the Preparation of Gram Negative Bacterial Genomic DNA:

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter

6.	Analytical balance
7.	Water distiller
8.	Vortex

The chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	Sodium hydroxide
2.	Glacial acetic acid
3.	SDS
4.	Di, sodium ethylene diamine tetra acetic acid (Na ₂ EDTA)
5.	Tris Base
6.	Chloroform
7.	Hydrochloric acid
8.	Ethanol
9.	Sodium chloride

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Glass centrifuge tubes
5.	Racks for Eppendorf tubes
6.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
7.	Biohazard Container
8.	Pipette

Preparation of Solution:

1. Lysis Buffer Solution:

It is prepared by:

1. Dissolve 0.242 g of Tris Base, 57 μL of Glacial acetic acid and 1 ml of 0.5 M EDTA solution.
2. To the above solution add 0.02 g of sodium acetate and 0.1 g of SDS.
3. The volume make up to 8 ml by distill water. Adjust the pH to (7.8), the volume complete to 10 ml by distill water.
4. Mix well, dissolve the solution with heating in water bath to assist dissolution.

2. 5 M NaCl Solution:

It is prepared by dissolving 2.92 g of Sodium Chloride in 10 ml of distill water.

3. Ethanol 70 %:

comprise of **70 ml** absolute ethanol and **30 ml** of distill water.

4. Ethanol 100 %:

Ready to use.

5. Chloroform:

Ready to use.

The Procedure:

1. 1.5 ml of a saturated culture was harvested with centrifugation for 5 min at 14,000 rpm.

2. The cell pellet was resuspended and lysed in **200 μL** of lysis buffer (40 mM Tris-acetate pH 7.8, 1 Mm EDTA, 20 mM sodium-acetate, 1 % SDS) by vigorous pipetting.
3. To remove most proteins and cell debris, **66 μL** of 5M NaCl solution was added and mixed well, and then the viscous mixture was centrifuged at 12,000 rpm for 10 min at 4°C.
4. After transferring the clear supernatant into a new Eppendorf tube, an equal volume of chloroform was added, Then mix on a rotating mixer for **15 min** when a milky solution was completely formed.
5. Following centrifugating at 14,000 rpm for 5 min, The supernatant is then removed to another Eppendorf tube and **double volume** of 100% EtOH was added.
6. The tubes were inverted **5-6 times** gently to **precipitated DNA**, then centrifuged at **10.000 rpm** for **3 minutes**.
7. The supernatant from this last step discarded and 1 ml of **ethanol (70%)** was added, again tubes centrifuged **10.000 rpm** for **3 minutes**.
8. **Repeating Step 7.**
9. The supernatant discarded and **100 μL** of distil water was added to re-suspended (dissolved) **DNA**.
10. The stock **DNA** is kept frozen at (-20°C) until use.

Genomic DNA Extraction Method From Gram Positive Bacterial

DNA extraction was according to Roeder and Broda (1987) and Thottappilly et al. (1999) with some modification.

The Instruments:

The following equipment's are used throughout this procedure:

No	Instrument
1.	Magnetic stirrer
2.	Shaking water bath
3.	Portable autoclave
4.	Microfuge
5.	pH-meter
6.	Analytical balance
7.	Vortex

The Chemicals:

The chemicals are used in this procedure:

No	Chemical
1.	CTAB
2.	2-mercaptoethanol
3.	Sodium chloride
4.	Di, sodium ethylene diamin tetra acetic acid (Na ₂ EDTA)
5.	Tris-Hcl
6.	SDS
7.	Hydrochloric acid
8.	Ethanol
9.	Sodium hydroxide

The Equipments:

The following instruments are used in this procedure:

No.	Equipment
1.	Rack for tubes
2.	Tips with different sizes
3.	Eppendorf tubes (1.5 size)
4.	Racks for Eppendorf tubes
5.	Micropipette (transfer pipette) (2-20µl, 10-100µl,100-1000µl)
6.	Biohazard Container
7.	Pipette

The Preparation of Buffers:

1. 2X CTAB buffer (50 mM Tris, pH 8.0; 0.7 mM NaCl; 10 mM EDTA; 2% CTAB; 0.1% 2-mercaptoethanol):

It is prepared by dissolving of **0.08 g Tris-HCL**, **0.001 g** of **NaCl**, **0.04 g** of **Na₂EDTA** **200 mg** of **CTAB**, and **10 µL** of 2-mercaptoethanol in **6 ml** of distill water. Adjust all to **pH (8)** using NaOH and make up volume to **10 ml** with distill water, sterilize by autoclaving.

2. 20% sodium dodecyl sulfate:

It is prepared by dissolving **2g** of **SDS** in **10ml** of distill water with heating at **68C** to assist dissolution.

3. Phenol : chloroform : isoamyl alcohol (25:24:1):

It is prepared by mixing 25 ml of saturated phenol (pH 8), 24 ml of chloroform and 1 ml of isoamyl alcohol in dark bottle.

4. Absolute ethanol:

Ready to use

5. 70% ethanol:

comprise of **7 ml** absolute ethanol and **3 ml** of distill water.

6. TE buffer:

It is prepared by adding 2.5ml of 1M Tris pH (8.0), to 0.5 ml of 0.5M EDTA pH (8.0), the volume made up to 250ml by distill water, autoclave and store at 4°C.

Procedure:

1. Transfer several colonies of gram positive bacterial isolates to 10 ml of nutrient broth and incubate for 24 hours at 37 °C.
2. Into the 1.5 ml eppendorf tube, pipette 1.5 ml of the overnight bacterial culture and spun at 14.000 rpm for 3 min in microcentrifuge. then washing with 200 µL TE buffer.
3. The washed bacterial cells suspension in 200 µL of 2xCTAB buffer, mixing briefly by vortex.
4. Add 100 µL of 20% SDS and incubated at 65°C for 20 min.
5. After incubation, add 200 µL of phenol : chloroform : isoamyl alcohol with gentle mixing, then centrifuged at 10000 rpm for 10 min.
6. Transfer the aqueous phase (upper layer) to new eppendorf tubes and add an equal volume of phenol : chloroform : isoamyl alcohol. The tubes mixing gently on a rotary mixer for 5 minutes, followed by centrifugation at 10000 rpm for 10 minutes.
7. Transfer the aqueous phase to another new eppendorf tubes and double volume of -20°C absolute ethanol, to precipitate DNA.
8. The tubes were inverted 5-6 times gently to precipitated DNA, then centrifuged at 10.000 rpm for 10 minutes. The supernatant from this last step discarded, then add 1 ml of ethanol (70%), again tubes centrifuged 10.000 rpm for 5 minutes.
9. The supernatant discarded and 100 µL of distil water, then add to re-suspended (dissolved) **DNA**.
10. The stock **DNA** is kept frozen at (-20°C) until use.

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Minipreparation or Rapid Plasmid Isolation Procedure:

1. Cultures from transformed colonies will be grown overnight in 2 mls of selective medium (LB + 50 ug/ml Ampicillin) in sterile 18 mm test tubes at 37°C and 250 rpms in the shaking incubator.
2. Pipette 1.5 mls of bacterial culture into eppendorf tubes and centrifuge cells in the Eppendorf centrifuge at maximum speed for 60 seconds to pellet cells.
3. Resuspend pellet thoroughly in 100 ul Lysozyme buffer, (Solution I). Note: Incomplete suspension of the cell culture will greatly reduce your plasmid yield.
4. Add a tiny bit of lysozyme on the end of a toothpick to the suspended cells and swirl toothpick in the tube to mix.
5. Incubate at RT° for 20 minutes.
6. Add 200 ul of lytic mix (Solution II) and MIX GENTLY by inversion.

Note: use freshly prepared lytic mix and treat the mixture gently from this point on. Rough handling will shear the chromosomal DNA into small pieces which will then contaminate your plasmid by co-precipitating with the plasmid.
7. Place on ice 5 minutes.
8. Add 150 ul of Potassium Acetate Buffer (Solution III) MIX GENTLY by inversion.
9. Place in ice water for 10 minutes to precipitate large cellular constituents such as cell wall, membrane proteins and chromosomal DNA.
10. Spin in the Eppendorf centrifuge 5 minutes.

11. Decant supernatant to clean eppendorf tube by pouring or by transferring with a Pasteur pipette.

12. Fill tube with RT° 95% EtOH-mix by inversion, and centrifuge 5 minutes.

13. Pour off ethanol, keep pellet, and dry pellet under vacuum in the Vacuum centrifuge evaporator for 10 minutes.

Note: It is important to remove all ethanol because it inhibits the activity of enzymes.

14. Add 30 ul of sterile nanopure water to the DNA pellet and mix well by vortexing for 1 minute to reconstitute the plasmid DNA.

15. Place the plasmid in a well -labelled tube at -20° C for long term storage.

Buffers for Alkaline Lysis DNA Extraction:

1. Solution I - 100 mls of Lysozyme buffer:

It is consist of:

- 50mM glucose (1.8 mls of 50% glucose).
- 25mM Tris HCl (2.5 mls of 1 M Tris, pH 7.6).
- 10 mM EDTA (2 mls of 0.5 M EDTA, pH 8).
- 94 mls of ddH₂O

2. Solution II - 10 mls of Lytic mix (make fresh each time):

It is consist of:

- 0.2 N NaOH final concentration (from 10 M Stock) 100 uL 10 M NaOH.
- 1% SDS (1ml 10% SDS).

- 8.9 mls ddH₂O.

3. Solution III - Potassium acetate buffer:

- 5 M K Acetate 60 ml.
- Glacial acetic acid 11.5 ml.
- H₂O 28.5 ml.

Note: Solution III will always be used cold, store at 4° C or on ice.

4. 5 M Potassium acetate (F.W. 98.14):

It is prepared by dissolving 49.07 g of Potassium acetate in 70 ml H₂O the volume is made up to 100ml by distill water Sterilize by autoclaving (15 minutes).

Glacial acetic acid:

Molecular Formula C₂H₄O₂

Molecular Weight 60.05

Isolation of plasmid DNA-small scale method:

World Journal of Microbiology & Biotechnology 19: 239–241, 2003.

Amjad B. Khalil^{1,*}, Ghandi H. Anfoka¹ and Salwa Bdour²

The described procedure is for 1.5 ml aliquots of an overnight culture of thermophilic bacteria or fresh single colonies from plates. Eppendorf-type tubes (1.5 ml) were used throughout and all centrifugations were carried out at maximum speed capable of generating 6–8000 rpm.

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The Preparation of Buffers:

1. Lysis Buffer: 10 ml

It is prepared by dissolving 0.09 g glucose, 0.035 g EDTA, 0.03 g Trizma base, 2.5 g sucrose) in 6 ml of distill water . Adjust the pH to (8.0). The volume is made up to 10 ml by the distill water and sterilize by autoclaving.

2. Lysozyme:

It is prepared by dissolving 1 mg of Lysozyme in 1 ml of distill water.

3. alkaline SDS Buffer: 10 ml (make fresh each time).

It is prepared by dissolving of 0.8 g NaOH and 1 g SDS in 10 ml distill water.

4. 7.5 M ammonium acetate:

It is prepared by dissolving 5.782 g of ammonium acetate in 10ml distilled water sterilize by filtration (micropore filter) with 0.20 μm .

5. Isopropanol: Ready to use.

6. TE buffer (10mM Tris–HCl, 1 mM EDTA, pH 8.0).

It is prepared by dissolving 0.016 g of Tris-HCl and 0.004 g of Na₂EDTA in 7 ml of D.W, Adjust the pH to (8.0). The volume is made up to 10 ml by the distill water and sterilize by autoclaving.

PROCEDURE:

1. The 1.5 ml culture aliquots were centrifuged for 2 min, the supernatant was decanted.
2. Pellet resuspended in 100 μL of lysis buffer.
3. Lysozyme was added at 1 mg/ml. The suspension was incubated at 37 °C for 15 min.
4. 100 μL of alkaline SDS was added to the mixture gently, and it was kept on ice for 15 min.

Note: Rough handling will shear the chromosomal DNA into small pieces which will then contaminate your plasmid by co-precipitating with the plasmid.

5. 100 μL of 7.5 M ammonium acetate was then added to the mixture and it was kept on ice for 15 min.
6. The mixture was centrifuged at 6000 rpm for 10 min. The supernatant was transferred to a new sterile centrifuge tube.
7. 400 μL of isopropanol was added to the new tube, mixed and kept at room temperature for 15 min.
8. Centrifuged at 7000 rpm for 10 min. The pellet was resuspended in 100 μL 2 M ammonium acetate and kept on ice for 15 min. The supernatant was saved and 200 μL of isopropanol was added to it, it was kept at room temperature for 10 min, and centrifuged at 8000 rpm for 15 min.
9. The pellet was dried for 15 min at room temperature, and then resuspended in 50 μL of TE buffer.

DNA Extraction from Paraffin-Embedded Tissues

Hongxin Fan and Margaret L. Gulley

1. Introduction:

In routine histopathology, most tissues are fixed in formalin and embedded in paraffin for long-term preservation. DNA can be extracted from these tissues for subsequent molecular analysis by amplification methods. We describe herein a protocol for DNA preparation from paraffin-embedded tissues based on published procedures (1–3). In brief, tissue sections are placed into microfuge tubes, then deparaffinized with xylene. The xylene is removed with ethanol washes, and the tissue is treated with proteinase K to make DNA available for amplification.

This protocol is simple, but there are several factors that influence the success of subsequent DNA amplification assays, including the type of fixative that is used, the duration of fixation, the age of the paraffin block, and the length of the DNA segment to be amplified (*see Note 1*). Ethanol fixation preserves DNA much better than does formalin. Formalin fixation randomly chops DNA in a duration-dependent manner, resulting in partial degradation. Even more severe DNA degradation occurs in bone samples subjected to acid decalcification. Because of this degradation, formalin-fixed tissue is not suitable for Southern blot analysis or for amplification of large DNA segments. Nevertheless, polymerase chain reaction (PCR) amplification of segments ranging up to 1300 bp has been reported (2), and consistent amplification of segments up to 300 bp is commonly achieved from archival fixed tissues.

Be aware that partial degradation of DNA may result in sampling bias, and therefore results should be interpreted with caution. For example, amplification of DNA from one cell may produce a PCR product that is not representative of the entire population of cells in the tissue. For this reason, it is wise to run tests in duplicate and always with appropriate controls.

2. Materials:

2.1. Reagents:

1. Xylene.
2. 100% Ethanol.
3. Proteinase K stock solution (20 mg/mL).
4. TEN buffer: 10 mM Tris-HCl, pH 8.0, 1 mM EDTA, pH 8.0, 20 mM NaCl.

5. 10% SDS.

2.2. Equipment:

1. Microcentrifuge.
2. Heating block or water bath to hold 1.5-mL Eppendorf tubes.
3. Microtome.

3. Methods:

Preparation of paraffin-embedded tissue for DNA amplification involves several manual manipulations; therefore precautions should be taken to avoid contamination, such as changing gloves frequently. When opening and closing microfuge tubes, do not touch the rim or inside of the cap. (Many laboratory scientists place a fresh gauze square over their thumb when opening a tube, or use a cap-opener device or screw-top lids.) Appropriate negative controls must be used to alert for contamination.

3.1. Cutting and Deparaffinizing Sections:

1. Use a microtome to cut five sections, 5–20 μm thick, from a paraffin block, and place these directly into a 1.5-mL microfuge tube. The thickness of the sections depends on the size of the biopsy. For a small biopsy (up to 3 mm), 20- μm thick sections may be required, whereas a large biopsy requires only 5- μm thick sections.

Although multiple thin sections can be placed in a single tube, fewer thick sections are more practical for processing. See **Note. 2** for special precautions against contamination.

2. Add 800 μL of xylene to each tube, close, mix by gentle vortexing, and then incubate at room temperature for 10 min. Pellet the tissue by centrifugation for 5 min in a

microfuge at full speed (20,000 rpm). Carefully remove and discard the supernatant using a pipet; do not disturb the tissue. If any translucent white paraffin remains, repeat the xylene wash one to two more times.

3. Add 800 μL of 100% ethanol to each tube, close the lid, and mix by inverting. Pellet the tissue by centrifugation for 5 min in a microfuge at full speed (20,000 rpm), and carefully remove and discard the supernatant with a pipet. Repeat the ethanol wash one more time, and remove as much supernatant as possible.
4. Open the tubes and let the residual ethanol evaporate by incubating in a dry heat block at 55°C for 15–30 min or until the sample is completely dry. (Speed vacuum drying is not recommended because of the risk of contamination.)

3.2. Proteinase K Digestion

12. To the dried tissue samples add 1 ml of TEN buffer containing 200 $\mu\text{g}/\text{mL}$ of proteinase K (prepared by mixing 10 μL of proteinase K stock solution in 1 ml of TEN buffer). and followed by **50 μL of 10% SDS**, mixed by vortexing briefly for 1 min until the mixture was re-suspended and well dissolved by gentle pipeting.
2. Close the tubes and incubate at 55°C for 2 hours. Gently invert the tubes a few times during incubation.
3. Spin briefly to remove any liquid from the cap. Cover the caps tightly with cap locks (PGC Scientifics, Gaithersburg, MD) to prevent them from popping open during high-temperature incubation. Incubate in a 95°C heat block for exactly 10 min to inactivate the proteinase K (*see Note 3*). Pellet the tissue in a microfuge at full speed

(20,000 rpm) for 10 min, and then transfer the supernatant to a clean tube and discard the pellet.

4. Add double volume from ice-cold ethanol and invert gently to allow DNA to precipitate. Centrifuge at 20,000 rpm for **10** minutes, discard the supernatant.
5. Wash the pellet by add **500 µL** of **70 %** ethanol. Centrifuge at **20,000** rpm for **10** minutes.
6. Resuspend in **50 µL** TE buffer.
7. Store DNA at -20°C , and avoid thawing and refreezing. (Freshly prepared samples are more efficiently amplified than those stored frozen, perhaps because the freeze-thaw cycle damages DNA.)

4. Notes:

1. The most important factor affecting DNA quality is the type of fixative employed and the duration of fixation. Tissues fixed between 12 and 24 h in ethanol, acetone, Omnifix, or 10% buffered formalin usually yield good-quality DNA; but B-5, Zenker's, or Bouin's solutions, or duration of fixation longer than 5 d, are poor prognostic factors for PCR productivity (2,4). Prior studies showed that when the length of α -globin amplification product increased (175, 324, and 676 bp), the percentage of fixed tissues containing amplifiable DNA decreased (100, 69, and 45%) (5). And when the age of a block increased, PCR productivity decreased (6,7), although some blocks stored for more than 40 yr were successfully studied (8).
2. It is important that no tissue be carried over from one case to the next during microtomy. Between each block, shift to a fresh part of the blade. Use smoothedge

rather than toothed forceps for transporting sections into tubes. Do not allow bleach to come in contact with the tissue or the DNA will be destroyed. Ice cubes used to cool a block should be discarded between cases.

3. Longer incubation at 95°C may damage DNA, whereas shorter incubation may not fully inactivate proteinase K. Additional time is needed for volumes >500 µL.

4. The optimal amount of template for an amplification reaction depends on numerous factors specific to each sample, such as DNA concentration and presence of inhibitors. It is useful to test several concentrations of each template (e.g., 1 and 10 µL of template per 100-µL PCR). Large tissues may necessitate the use of a smaller fraction of the template (e.g., 0.1 µL). The amount of template that is “tolerated” in a PCR may be affected by residual fixation chemicals or paraffin, excessive tissue debris, and other factors. If the first attempt fails, a 10-fold dilution will often reduce inhibitors while still retaining enough DNA to allow amplification.

5. Amplification of DNA prepared from paraffin-embedded tissue is less efficient than amplification of DNA from fresh or frozen tissues. To compensate for this reduced efficiency, consider modifying the thermocycling parameters by increasing the number of cycles and lengthening the duration at each temperature within the cycle (1).

Reference :

Methods in Molecular Medicine, vol. 49 : Molecular Pathology Protocols

Edited by: A. A. Killeen © Humana Press Inc., Totowa, NJ

DNA Extraction From Plasma and Serum:

1. The chemicals:

No.	Chemicals
1.	Di, sodium ethylene diamin tetra acetic acid (Na_2EDTA)
2.	Tris-base
3.	Tris-Hcl
4.	Protienase K
5.	Hydrochloric acid
6.	SDS
7.	Ethanol
8.	Glycogen
9	Sodium hydroxide
10.	Ammonium acetate
11.	Phenol
12.	Chloroform

2. The Equipments:

No.	Equipments
1.	Magnetic stirrer
2.	Portable autoclave
3.	Centrifuge
4.	Microfuge
5.	Shaking water bath
6.	pH-meter
7.	Analytical balance
8.	Water distiller
9	Vortex
10.	MicroCentrifuge

3. The Instruments:

No.	Instrument
1.	Micropipette (transfer pipette) (2-20 μ l, 10-100 μ l, 100-1000 μ L).
2.	Biohazard Container
3.	Biohazard Bags
4.	Glass centrifuge tubes
5.	Pasteur pipettes
6.	Eppendorf tubes
7.	Tips with different size
8.	Pipettes (1ml, 5ml, 10ml)

Preparation of Stock Solutions:

1. 1M Tris HCl pH 8.0 :

Prepare by dissolving 60.55g of Tris base in 400 ml of distill water . Concentrate HCl are added to adjust the pH to (8.0). The volume is made up to 500ml by the distill water and sterilize by autoclaving.

2. 0.5 M EDTA pH 8.0 :

Prepare by dissolving 93.05g of EDTA in 400ml of distill water and stir vigorously on a magnetic stirrer . The pH is adjust to (8.0) using (10N) NaOH. The volume is made up to 500ml by distill water and sterilize by autoclaving.

3. 10M NaOH :

It is prepared by dissolving 40g of NaOH in 80ml of distill water the volume is made up to 100ml by distill water.

4. Phenol :

It is prepared by melting crystalline phenol at 68C water bath and then overlaid with an equal volume of 1M Tris pH (8.0) and is allowed to stand overnight at 4C. The 1M Tris is then discard to be replace by 0.1M Tris pH (8.0), store at 4C for up to 4 weeks .

5. 10% SDS Solution :

It is prepared by dissolving 100g of SDS in 900ml of distill water with heating at 68°C to assist dissolution adjust volume to one liter.

6. 7.5M Ammonium acetate :

It is prepared by dissolving 5.782g Ammonium acetate in 10ml distilled water sterilize by filtration (micropore filter) with 0.20 μm .

7. Absolute ethanol (100%):

Ready to use

8. Ethanol 70% :

Comprise of 70 ml absolute ethanol and 30 ml of distill water.

9. Glycogen (10 mg/mL):

It is prepared by dissolving 10g of Glycogen in 1ml of distill water store at 4°C until use.

10. Proteinase K 5 mg/mL:

It is prepared by dissolving 5 mg of Proteinase K in 1ml of distill water store at -20 °C until use.

Preparation of Buffers (Work) Solutions:

1. TE Buffer :

It is prepared by adding 2.5ml of 1M Tris pH (8.0), to 0.5 ml of 0.5M EDTA pH (8.0) , the volume made up to 250ml by distill water, autoclave and store at 4°C.

2. Phenol : chloroform (1:1, v/v):

It is prepared by mixing one volume of equilibrated phenol pH (8.0) with one volume of chloroform This solution store in a light-tight bottle at 4°C.

3. 10X SDS/Proteinase K:

It is consist of 10 g/100 mL lauryl sulfate (SDS), 5 mg/mL proteinase K.

For each sample of serum or plasma 1.5 mL of 1X SDS/Proteinase K solution is prepared by mixing the following volumes:

1. 150 μL of SDS (10%), (As final concentration 1%).

2. 7.5 μL of Proteinase K (5 mg/mL), (As final concentration 50 $\mu\text{g/mL}$).

3. The volume is made up to 1.5 ml (Final volume) by add 1342.5 of distill water.

Procedure:

1. Place 1.5 mL serum or plasma into a 10-mL centrifuge tube.
2. Add 1.5 mL of 1X SDS/proteinase K solution in the tube containing the serum and mix well.
3. Digest overnight at 55°C in water bath.
4. Add 3 mL phenol :chloroform solution.
5. Vortex for 30 s and centrifuge for 10 min at 1000g using a swing-out rotor.
6. Transfer aqueous layer to fresh tube and repeat **steps 4 and 5**.
7. Transfer aqueous layer to fresh tube and add 5 µL glycogen (10 mg/mL), 1 mL of 7.5M ammonium acetate, and 8 mL of 100% ethanol.
8. Mix by inverting and centrifuge at 2500g for 40 min.
9. Carefully remove supernatant and wash pellet in 10 mL of 70% ethanol.
10. Centrifuge at 2500g for 10 min. Carefully remove last traces of ethanol and allow to air-dry for 10 min before redissolving the pellet in 100 µL TE.