




Enterobacteriaceae



Assist. Prof. Dr. Ihsan Muneer Ahmed
Department of Environmental Health

27.2.2025

1

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Enterobacteriaceae



The most important species:

- 1- *Escherichia coli*
- 2- *Klebsiella* spp.
- 3- *Salmonella* spp.
- 4- *Proteus* spp.



2

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Enterobacteriaceae



Enterobacteriaceae general characteristics:

- 1- Gram-negative non-spore forming bacteria.
- 2- Medium-sized (0.4-0.6×2-3µm).
- 3- The motile species have peritrichous arrangement of flagella.
- 4- Facultative anaerobic.
- 5- Catalase - positive and oxidase - negative.
- 6- Reduce nitrate to nitrite.
- 7- Able to grow on non-enriched media such as nutrient agar.

3



Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Enterobacteriaceae on MacConkey agar:

Lactose fermentation in MacConkey agar:

- 1- The colonies of lactose fermenters are **pink** due to acid production from lactose (*E. coli*, *Klebsiella*).
- 2- The colonies of non-lactose fermenters have **pale** and are alkaline due to utilization of peptones in the medium (*Salmonella*, *Shigella*, *Proteus*).



4

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Enterobacteriaceae classification according to biochemical reaction tests (Table 1):

I	Indole production
M	Methyl red
V	Voges – Proskauer
C	Citrate utilization



IMVC

H	H ₂ S production
U	Urease
M	Motility



HUM

5

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Table 1: Biochemical reactions for main members of the Enterobacteriaceae

Genus	I	M	V	C	H	U	M
<i>E. coli</i>	+	+	-	-	-	-	+
<i>Klebsiella</i>	-	-	+	+	-	+	-
<i>Salmonella</i>	-	+	-	+	+	-	+
<i>Proteus</i>	V	+	-	-	+	+	+

v = variable

6

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Diagnosis of Enterobacteriaceae



Specimens for bacterial examination should be cultured aerobically on MacConkey agar at 37°C for 24 hours.

General identification criteria for isolates include:

- 1- Gram-negative rods.
- 2- Catalase-positive, oxidase-negative.
- 3- Growth and appearance on MacConkey agar (lactose fermenter or non-lactose fermenter).
- 4- Biochemical profile including IMVC and HUM)for definitive identification.

7

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Diseases produced by Enterobacteriaceae



ORGANISM	DISEASES
Lactose Positive	
<i>Escherichia coli</i>	Diarrhea, sepsis, urinary tract infection, neonatal meningitis
<i>Klebsiella, Serratia, Citrobacter, Enterobacter</i> species	Opportunistic infections (e.g., pneumonia, sepsis, neonatal meningitis)
Lactose Negative	
<i>Salmonella</i> species	Diarrhea, typhoid fever, bacteremia; localized infections in bone, meninges, liver
<i>Shigella</i> species	Dysentery
<i>Proteus</i> species	Urinary tract infection
<i>Yersinia</i> species	Plague, diarrhea, mesenteric lymphadenitis

8

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



1- *Escherichia coli*




Morphology and Staining

Cultural Characteristics

Biochemical Tests

Diagnosis



9

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025




Morphology and Staining



Escherichia coli

Theodor Escherich, 1885.



Gram-negative, medium-sized rod (3-1 × 0.6 μm).

10


Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Cultural Characteristics



Escherichia coli on MacConkey agar **ferment lactose** and produce **pink** colonies.



11

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



E. coli on Eosin Methylene Blue agar (EMB) produce the characteristic phenomena greenish color called metallic sheen.



Metallic sheen

12

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Biochemical tests



Species	Indole production	Methyl red	Voges – Proskauer	Citrate	Hydrogen sulphide	Urease	Motility (36 °C)	Gelatin liquefaction	Lactose	Maltose	Mannitol	Sucrose
<i>E. coli</i>	+	+	-	-	-	-	(+)	-	+	+	+	d

+ = 90 – 100% strain positive
(+) = 76 – 89% Positive

- = 0 – 10% Positive
d = 26 – 75% Positive

13

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Escherichia coli



+

+

-

-

+

-

A/A/-



I

M

V

C

N

U

H

I = Indole, M = Methyl red reaction, V = Voges – Proskauer reaction,
C = Citrate utilization, N = Nitrate reduction, U = Urease,
H = Hydrogen sulphide production (H₂S) test, A = Acid,

14

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025




Diagnosis



Identification criteria for *E. coli*.

- On MacConkey agar, colonies are bright pink.
- Using IMVIC tests for confirmation (Table.1).
- On EMB agar, the colonies of *E. coli* have a metallic sheen.



15

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



2- *Klebsiella* spp.




Morphology and Staining

Cultural Characteristics

Biochemical Tests

Diagnosis



16

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025





Morphology and Staining



Klebsiella pneumoniae

1-Gram -negative, medium-sized rods $2-1 \times 0.5-0.8 \mu\text{m}$.

2-Non-motile (because of capsule), non-spore forming bacteria.



Capsule stain

17


Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Cultural Characteristics



Klebsiella pneumoniae on MacConkey agar is a lactose-fermenter with characteristic large mucoid pink colonies.



18

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Biochemical tests



Species	Indole production	Methyl red	Voges – Proskauer	Citrate	Hydrogen sulphide	Urease	Motility (36 °C)	Gelatin liquefaction	Lactose	Maltose	Mannitol	Sucrose
<i>K. pneumoniae</i>	-	(-)	+	+	-	+	-	-	+	+	+	+

+ = 90 – 100 % strain positive

- = 0 – 10 % Positive

(-) = 11 – 25 % Positive

19

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Klebsiella pneumoniae



—

—

+

+

+

+

K/A/-



I

M

V

C

N

U

H

I = Indole, M = Methyl red reaction, V = Voges – Proskauer reaction, C = Citrate utilization, N = Nitrate reduction, U = Urease, H = Hydrogen sulfide production (H₂S) test, K= Alkaline, A = Acid,

20

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025





Diagnosis



Identification criteria for *Klebsiella pneumoniae*.

- On MacConkey agar, colonies are mucoïd pink.
- Using IMVIC tests for confirmation (urease +ve) (Table.1)
- Using capsule stain (negative stain), *Klebsiella pneumoniae* have clear capsule.



21

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



3- *Salmonella* spp.




Species

Morphology and Staining

Cultural Characteristics

Biochemical Tests

Diagnosis



22

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Species



All *Salmonella* species are divided into 2 main species only:

1. *Salmonella enterica*

- *Salmonella enterica* subsp. *enterica* ser. *Typhi*
(*Salmonella typhi*)
- *Salmonella enterica* subsp. *enterica* ser. *Typhimurium*
(*Salmonella typhimurium*)

2. *Salmonella bongori*

23

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025




Morphology and Staining



Salmonella spp.

Daniel Elmer Salmon, an American veterinary surgeon, 1900.



Gram-negative, medium-sized rod (3-1 × 0.6 μm).

24



Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Cultural Characteristics




1. On MacConkey agar *Salmonella* spp. produce pale non-lactose fermenter colonies.
2. On Salmonella Shigella agar is a *Salmonella* spp. produce characteristic black center colonies.
3. Tetrathionate broth and selenite broth enrichment broth for salmonella growth.




25

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



MacConkey Agar
Pale colonies



Salmonella Shigella Agar
Black center colonies

26

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Biochemical tests



Species	Indole production	Methyl red	Voges – Proskauer	Citrate	Hydrogen sulphide	Urease	Motility (36 °C)	Gelatin liquefaction	Lactose	Maltose	Mannitol	Sucrose
<i>Salmonella</i>	-	+	-	+	+	-	+	-	-	+	+	-

+ = 90 – 100 % strain positive

- = 0 – 10 % Positive

27

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Salmonella



—

+

—

+

+

—

K/A/H₂S



I

M

V

C

N

U

H

I = Indole, M = Methyl red reaction, V = Voges – Proskauer reaction, C = Citrate utilization, N = Nitrate reduction, U = Urease, H = Hydrogen sulfide production (H₂S) test, K = Alkaline, A = Acid,

28

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025




Diagnosis



- Identification criteria for *Salmonella*.

- Specimens should be cultured in tetrathionate broth for preliminary enrichment and subculture in selective plating media such as (Salmonella-Shigella agar, Brilliant green agar).
- On MacConkey agar, colonies are pale.
- On Salmonella Shigella agar produce black center colonies.
- Using IMVIC tests for confirmation especially, (TSI +ve) (Table.1).



29

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



4- *Proteus* spp.



Morphology and Staining

Cultural Characteristics

Biochemical Tests

Diagnosis



30

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Morphology and Staining





Proteus mirabilis

Proteus vulgaris

1-Gram - negative, small rods $1-3 \times 0.5 \mu\text{m}$

2-Actively motile peritrichous flagella, non capsulate.



31


Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Cultural Characteristics



- 1- On MacConkey agar, *Proteus* do not ferment lactose and form pale colonies.
- 2- On blood agar, *Proteus* spp. produce swarming phenomena due to vigorous motility, it can be inhibited by increasing the concentration of agar in the medium.
- 3- On brain heart infusion agar, *Proteus* spp. growth (swarming phenomenon).



MacConkey agar



Blood agar



BHI agar

32


Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Swarming phenomena for *Proteus*



Blood agar



BHI agar

33

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Biochemical tests




Species	Indole production	Methyl red	Voges - proskauer	Citrate	Hydrogen sulphide	Urease	Motility (36 °C)	Gelatin liquefaction	Lactose	Maltose	Mannitol	Sucrose
<i>Pr. mirabilis</i>	-	+	(-)	d	+	+	+	+	-	-	-	(-)
<i>Pr. vulgaris</i>	+	+	-	(-)	+	+	+	+	-	+	-	+

+ = 90 – 100% strain positive
(-) = 11 – 25 % Positive


- = 0 – 10% Positive
d = 26 – 75% Positive


34

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025





Proteus





K/A/H₂S





I

M

V

C

N


U

H


I = Indole, **M** = Methyl red reaction, **V** =Voges – Proskauer reaction,
C = Citrate utilization, **N** = Nitrate reduction, **U** = Urease,
H = Hydrogen sulfide production (H₂S) test, **K**= Alkaline, **A** = Acid,

35

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



Diagnosis



Identification criteria for *Proteus* spp.

- On MacConkey agar, colonies are pale.
- Using IMVIC tests for confirmation (Table.1). (urease +ve)
- On BHI agar and blood agar, *Proteus* produce swarming phenomena.

36

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025



37

Practical Medical Bacteriology – 2nd stage - Department of Environmental Health © 2024-2025